



V

INFINI<sup>®</sup>  
MEDICAL DEVICES



**V<sup>E</sup>**

**An innovative dermal filling and conditioning product designed for the delicate skin around the eye area.**

**5ml**

V<sup>E</sup> is an advanced tissue stimulator that combines the benefits of fragmented hyaluronic acid and an amino acid complex to provide deep revitalization and structural support to the delicate skin around the eye area. The formula is designed to work three-dimensionally to hydrate, delay the aging process and eliminate bruising, dark circles and sagging of the eyelid area.

A key element of this product is the presence of hyaluronic acid fragments (HAF) with a low molecular weight (20-38 monomers), which allows for effective stimulation of CD-44 receptors in fibroblasts. In addition, the amino acids it contains intensify the process of skin protein remodeling. As a result, the product actively supports the processes of collagen, elastin and glycosaminoglycan synthesis. The result is a visible improvement of the eye contour along with intensive regeneration and smoothing of the skin.

**INGREDIENTS:**

**Fragmented hyaluronic acid:**  
**20-38 monomers.**

**Amino acid complex:**  
**Isoleucine, Leucine, Lysine, Proline, Valine,  
Glycine, Serine, Alanine, Cysteine**

**CE**  
**0373**

**PREMIERE:**  
**MAY 2024**



# V1

An advanced tissue stimulator for restoring skin tone and elasticity

5ml

An advanced tissue stimulator for restoring skin tone and elasticity. Its unique formula is based on a combination of fragmented hyaluronic acid, with lifting choline and amino acid complex. Choline, as a precursor of acetylcholine - an important neurotransmitter - plays a key role in maintaining proper muscle and skin tone. Its effect on acetylcholine synthesis contributes to the activation of cholinergic skin functions, which helps maintain a taut, well-defined facial contour and skin texture.

Hyaluronic acid fragments (HAF), with chain lengths ranging from 20 to 38 monomers, effectively stimulates CD-44 receptors in fibroblasts.

The effect is to intensify skin remodeling processes and increase the production of collagen, elastin, and glycosaminoglycans.

As a complement to the composition, an amino acid complex accelerates tissue regeneration and promotes the formation of new protein structures.

**INGREDIENTS:**

**Fragmented hyaluronic acid:**  
20-38 monomers

**Amino acid complex:**  
Serine, Alanine, Cysteine, Lysine, Proline, Glycine

**Choline:**  
Acetylcholine precursor



**PREMIERE:**  
MAY 2024



# V2

An innovative cellular stimulator that delivers revolutionary results in beautifying, skin rejuvenation and skin antioxidant protection.

5ml

The unique formula combines a powerful amino acid-vitamin complex, enriched with fragmented hyaluronic acid and potent antioxidant glutathione.

The key ingredient is fragmentary hyaluronic acid (20-38 monomers), which actively stimulates CD-44 receptors in fibroblasts, which promotes collagen synthesis and remodeling. This advanced formula not only rebuilds the skin's supporting structure, but also has antioxidant and firming effects, and guarantees long-lasting hydration and brightening of the skin.

**INGREDIENTS:**

**Fragmented hyaluronic acid:**  
20-38 monomers

**Amino acid collagen promoters:**  
Serine, Alanine, Cysteine, Leucine, Lysine, Proline, Valine, Glycine

**Reduced Glutathione**

**Stable form of vitamin:**  
Sodium Ascorbic Phosphate Buffer



**PREMIERE:**  
MAY 2024





# V5


SKIN STIMULATOR (H-HA) + (L-HA) + TREHALOSE 70mg / 2ml

A unique formula that combines hyaluronic acid with a special range of molecular weights and trehalose, an extremely effective stabilizing ingredient and creates a more intense moisturizing effect and greater, more powerful skin regeneration.

Hyaluronic acid brings instant hydration and volume, improving skin texture and geometry. At the same time, trehalose protects the acid molecules from enzymatic degradation, prolonging the skin's hydration effect.

Thanks to its anti-glycation properties, it helps to protect skin cells from oxidative stress, optimizing tissue regeneration and delaying the aging process.

The combined complex helps to improve skin condition and maintain its youthful appearance by hydrating, stimulating collagen and protecting against aging factors and environmental damage.



70mg / 2ml

CE  
1984

PREMIERE:  
MAY 2024



# V10

**Volumetric Stimulator**  
**365mg Poly-L-Lactic Acid**

V-10 Volumetric Stimulator uses powerful properties of Poly-L-lactic acid, specially designed to effectively lift the skin and correct a variety of skin defects such as wrinkles, furrows, scars and fatty tissue atrophy. Its unique formula is designed not only to even out the surface of the skin, but also to significantly improve its quality and density.

The action of Poly-L-lactic acid is focused on stimulating fibroblasts to produce collagen and replenish the atrophy of fatty tissue. This controlled biosynthesis of new collagen fibers leads to an increase in tone and volume of soft facial tissues, providing visible effects of rejuvenation and improvement of skin structure.

In terms of its action, it not only remodels facial contours, but also effectively promotes skin regeneration, increasing its thickness and elasticity. This dual action guarantees significant skin rejuvenation and restoration of lost facial volume and restoration of a harmonious appearance.


V10- is a product filled with freeze dried Poly-L-Lactic acid (PLLA) powder in a vial. PLLA is injected into the deep dermis or subcutaneous tissue to temporarily correct moderate to deep nasolabial fold wrinkles in adults through physical restoration.

**Vial contains:**  
**150mg of Poly-L-lactic acid**  
**90mg of Sodium Carboxymethylcellulose**  
**125mg of Mannitol**

**CE**  
**1434**

**PREMIERE:**  
**APRIL 2024**





# V20

**Advanced skin biovitalizer based on high-concentrated hyaluronic acid 40mg/2ml**

V 20 is an exclusive composition containing a uniquely formulated polymer of uncrosslinked hyaluronic acid at a high concentration of 40mg/2ml. It is characterized by a medium molecular weight in the range of 1,200 - 1,500 kDa and a unique spatial structure. Thanks to the advanced technologies implemented, it prevents the effects of the main factors of skin aging through its intensive stimulation and regeneration.

Thanks to its unique structure, it reconstructs the extracellular matrix, which provides an ideal environment for the biosynthesis of new collagen and elastin fibers.

The product has the ability to intensively bind water, which translates into long-term hydration of the skin. In addition, it promotes the processes of biological regeneration, which contributes to increased skin hydration, tone and firmness, while slowing down the signs of aging.



**PREMIERE:**  
MAY 2024



# V BODY

An intensive cellular stimulator enriched with fragmented hyaluronic acid and the amino acids lysine, glycine and proline designed specifically for body and facial firming treatments.

9ml

Fragmented hyaluronic acid, which is a polysaccharide with a molecular weight of 20-38 monomers, is characterized by its unique ability to retain water, which is key to maintaining optimal cell hydration and modulating the extracellular matrix.

Its unique feature also lies in the activation of CD-44 receptors, which promotes the biosynthesis of collagen and elastin, thereby promoting the structural reconstruction of the skin.

The presence of amino acids, fundamental components of proteins, further enhances regenerative processes, improving hydration, elasticity and overall skin condition. It is recommended for aesthetic treatments targeting sagging skin on the face and body.

Its use is indicated for the face, neck, décolleté, as well as other parts of the body such as the buttocks, hips, abdomen, inner side of the arms and thighs, where the skin requires intensive firming therapy.

**INGREDIENTS:**

**Fragmented hyaluronic acid:**  
20-38 monomers

**Amino acid collagen promoters:**  
Lysine, Glycine, Proline

CE  
0373

**PREMIERE:**  
MAY 2024





INFINI LAB S.R.L.  
Via Francesco Sforza n.15  
20122 Milano  
Italia

[infinifiller.com](https://www.infinifiller.com)

[contact@infinifiller.com](mailto:contact@infinifiller.com)

+39 028 974 58 80

REV 001/2024